# CANULE UROLOGIQUE

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### DESCRIPTION

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## INSTRUCTIONS FOR USE

### DESCRIPTION

These instructions for use must be read before using the urological cannula.  
Urological cannula is used for administration of medicated paste or cream preparations.  
Urological cannula is manufactured by **UNION PLASTIC** by injection of plastic materials compliant with food contact and/ or European Pharmacopeia. Urological cannulas are decontaminated with Ethylene Oxide.  
Some references may be colored by the use of a coloring master batch food contact compliant.  
Cannula is individually packaged in a sealed package.  
Cannula must be used with the drug for which it is prescribed and according to the pharmaceutical manufacturer's instructions for use.  
Cannula is a single use device and must be used for a single patient.

### USE

Cannula must be used in a healthy environment.  
Cannula must only be used by an health professional and in good hygienic conditions.  
Before use, it is necessary to wash your hand.

- 1. Preparation of drug's tube:**
  - Get the tube out its packaging,
  - Disinfect the tube and the cap with alcohol,
  - Unscrew the cap,
  - Turn the cap over and drill the top of the tube,
  - Pull the cap out.
- 2. Preparation of cannula's use:**
  - Wear sterile chirurgica gloves,
  - Open the cannula's packaging,
  - Screw the cannula on the drug's tube.
- 3. Administration to the patient:**
  - Can only be used by urological way according to pharmaceutical manufacturer's instructions for use.
  - Press the tube to dispense the drug.
- 4. Throw away the cannula:**
  - At the end of the handling, throw away the cannula and the tube according to the procedure in force concerning the disposal of infectious risk wastes.

### INDICATIONS FOR USE

Cannula must be used according to pharmaceutical manufacturer's instructions for use and to drug.  
It is a single use Medical Device.  
These instructions for use should be translated in drug's instructions for use in using country official language.  
Any serious incident related to the device must be notified to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

### CONTRAINDICATION

- Do not use cannula if its packaging is damaged,
- Do not use cannula if it is damaged,
- Do not use cannula for another application or with another drug than it is prescribed for,
- Do not use cannula for administration in another orifice,
- Do not use cannula when it is expired.
- Do not use cannula for a child weighing less than 2.4 kg.
- Do not use the cannula beyond 5 minutes of contact time with anal orifice.

### MEDICAL DEVICE ACCURACY

There is no claim accuracy for cannula.

### PRECAUTIONS OF USE

**DO NOT USE CANNULA IF IT IS DAMAGED OR EXPIRED.**  
Cannula is ready to be used, do not wash or sterilized it (do not sanitize or boil it).  
Cannula must be used by people who have no insanity of mind or impaired mental faculties and can read and understand the instructions for use.  
Cannula must only be used by a health professional.  
Check if cannula's packaging is not damaged.

### CONDITIONS OF STORAGE

Cannula must be kept in its own origin packaging at a temperature between 5°C and 50°C.  
After use, cannula and the drug's tube must be thrown away in appropriate way.

### MANUFACTURER INFORMATIONS

Expiry date on cannula's packaging.  
Claimed expiry date: 5 years.  
First placing date on the market: 2022.  
Last revision of the instructions for use: Revision D 2022/07/08  
Class 1 Medical Device.

Medical Device

Catalogue Reference

Batch Number

Unique Device Identification

Expiry date

Manufacturer

Do not reuse

Do not use if the packaging is damaged

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UNION PLASTIC SAS